

Using the METRIC Client

Robert Rommel

BioMedware

June 23, 2017



BioMedware
Geospatial Research and Software

Caveat

- a. We are very early in the development cycle
- b. Features are very limited
- c. Some bugs may exist

Navigating Through the Views

- Icons at the top handle main navigation
- You can move in any direction
- Hover over an icon will tell you what the view is



BioMedware
Geospatial Research and Software

The Project Tree

- ➔ Located in the top right
- ➔ Displays datasets that you have loaded into your session
- ➔ Displays results that you have calculated
- ➔ Right now, you can interact with these objects by right-clicking
- ➔ In the future you will be able to smartly drag/drop these objects and the display will be more informative

Using the Visualizations

- ➔ The displays are created automatically
- ➔ Map allows for layers, but currently only implement a single coloring scheme (analysis results are handled differently)
- ➔ Currently if you want to change a visualization, re-add the data to the visualization
- ➔ In the future you will be able to customize the appearance of the displays and change the data from the visualization itself

Running an Analysis

- ➔ Only appropriate data can be used (we aren't yet checking for variable types)
- ➔ Analyses are fairly swift, results are visualized automatically